



Connect Your iTero® to BioDent's Workflow

Save time and money when you send scans directly from your iTero® digital impression scanner. **In addition to all-ceramic restorations, BioDent accepts digital impressions for PFM and full-cast restorations.**

- Call iTero tech support at 800-577-8767, select option 2 and ask them to add BioDent as a lab partner.
- Select BioDent as your lab partner and submit your digital file.

If you have any questions, contact BioDent at 800-517-5250 and ask for digital impression technical support.

Looking for more information on our products?

A complete listing of restorative options and material descriptions, preparation and cementation guides and billing codes, visit biodentlaboratory.com/for-dentists

Save money on every digital submission
by choosing model free when possible.

SAVE \$20 On Your First Digital Case

Digital submissions are simple when you choose BioDent. Reference coupon code **GO_DIGITAL** on your Rx when you upload your first case and save \$20 on your first digital submission.

Coupon Code: GO_DIGITAL